## PHYSICIAN CONNECT All Providers

## Nuance PowerMic Mobile (Android)

Before you begin, make sure you have your Google account information ready.

1. On your Android device, go to the Google Play store.



2. Tap the **search** button (**magnifying glass**) at the top, right-hand side of your **Google Play** app, and search for the **Nuance PowerMic Mobile** app. Install the app.



3. On your Android device, open the browser (e.g., Google Chrome) and type this URL:

https://it.houstonmethodist.org/dragon/. The following Dragon resources/training page opens.



4. Scroll to the bottom of this webpage to **Quick Links**. Click on the **PowerMic Mobile Configuration** link for **Android** to complete the setup.



- 5. Once you click the link, tap **Open with**, **PowerMic Mobile**.
  - a. On the Default app selected pop-up window, select OK.
  - b. Accept the License Agreement.
  - c. Tap Start.



- 6. Enter your HM network ID/Epic user ID.
  - a. Tap Log In.
  - b. Select Allow to Allow PowerMic mobile to record audio.

10	4 22 6		94	<u>8</u> 13:35	
	Nuance' Powe Mobi	erMic le			
hr	mstjinxg				
		Log In			
3 Þ		1914 🗎 1429			
IMSTJN06	вникст		8	9 0	
	✓		1	o p	
			k	1	
D.		T+∎	m	a	
Allow PowerMic Mobile to record audio?					
	DENY	ALLOW			
(	٢				
×		H			

7. The **PowerMic Mobile** app will begin searching for a workstation to pair with. Make sure you are logged into Epic and have launched Dragon from the **Dragon Login** button on the **Epic toolbar**.



8. Make sure the icon on the **DragonBar** is a **cell phone**. If it's a **red microphone**, you must click the **DragonBar Menu** (flame) and hit Log Off.



9. In a few seconds, a window or splash screen pops up that allows you to change the mic source to **Nuance PowerMic Mobile**.

Nuance <sup>®</sup> Dragon Medical One		Desktop Application
	User:	
	Microphone	
	Microphone (PowerMicIII-NS)	
	Microphone (PowerMicIII-NS)	
	Nuance PowerMic Mobile	
	Specialty:	
	General Medicine	•
	ОК	Exit
		NUANCE
© 2011-2017 Nuance Communications, Inc. All rights reserved.		3.5.101.2632

10. Once you have logged into Dragon on your computer using the mic source as **Nuance PowerMic Mobile**, the app on your phone will pair with the computer and the **orange cell phone** icon on your **DragonBar** will turn **red**.



11. The app opens on your phone. Tap the **mic** button on the app to begin dictating. Once you tap the **mic** button on your app, the **red cell phone** icon on your computer will turn **green**. Start dictating with the cursor in the field/note section where you want the text to appear.

