

Caregility – vICU Login



Audience: vMD and vRN

Getting Started

LOGGING INTO CAREGILITY

- To login via desktop, go to <https://landing.caregility.net/>
- Caregility can be launched using Chrome and Safari only.
- If you're off the HM network and logging in, you'll be prompted to go through SSO.
- Enter Email.

Caregility Login

Log into your Account

Email

Continue

If you're having trouble accessing the system, please contact your company's administrator or contact:

Caregility Technical Support
Toll free: 866-300-0884
or email: techsupport@caregility.com

- HM username and password.

HOUSTON
Methodist
LEADING MEDICINE

Caregility SSO External

Please select the 'Public Computer' option if this is not a machine you use regularly, then enter your User ID below and click 'Submit' to access the system.


Username:

This is a public computer
 This is a private computer

Submit

Caregility – vICU Login

- Select the **Program** after logging in.



Caregility Login

Select Your Program


Program

HM vICU

Testing Mode

Select

- Select the **iConsult** icon from the home screen. If prompted, select **Allow** for camera and audio.




mobile.uheplatform.com wants to

Use your microphone

Use your camera

Allow Block

CALLING VICU WALL/CART



Clinician Guest Mobile Device **Wall/Cart** SIP/H.323

Bedside Favorites

HM Houston Methodist

HMCC Houston Methodist Continuing Care

HMCC LTACH 2 WEST ICU

HMCC LTACH 2 West ICU 228

Call in Privacy Mode

CALL

Caregility – vICU Login



- Call is initiated.



	Camera Control
	Camera Home
	Settings
	Bookmarks
	Self View
	Patient on Hold
	Microphone Mute
	Camera Mute
	Door Bell



HOW TO CALL A CLINICIAN – BRIDGE CONSULT

- Click **Participants**.
- Choose the **Invite Type**.

