

EPOprostenol (Veletri) (RESTRICTED) [1205]

Nursing

Nursing

- | | |
|--|--|
| [X] Notify ACCREDO for initiation of EPOprostenol therapy, must verify enrollment. Dosing weight may be different from current weight . For initiation of therapy, must verify enrollment with: ACCREDO 1-866-344-4874 or CAREMARK 1-877-356-5264 for financial responsibility for home therapy. | Routine, Until discontinued, Starting S |
| [X] EPOprostenol (Veletri) instructions for Nursing: | Routine, Until discontinued, Starting S
1. Do NOT infuse any other IV fluids or medications through Veletri line.
2. Do NOT draw blood from a lumen adjacent to Veletri line.
3. Do NOT STOP Veletri for ANY REASON without first contacting MD who ordered Veletri
4. DO NOT STOP Veletri for CODE without first contacting MD who ordered Veletri. |

IV Fluids

Medications

Medications (Single Response)

- | | |
|---|--|
| () EPOprostenol (VeLETri) straight dose | intravenous, for 24 Hours, daily at 1400 (time critical)
RESTRICTED to ICU Intensivists, Pulmonology and Cardiology specialists. Are you an ICU Intensivists, Pulmonology or Cardiology specialist or ordering on behalf of one?
Infusion Pump Rate (mL per 24 hours): |
|---|--|

Additional Orders