

Location: \_\_\_\_\_

**Nursing**

**IV Fluids**

**IV Fluids**

- ☐ **heparin 1000 unit/500 mL 0.9% sodium chloride** Dose: 100 mL Route: intravenous Frequency: once in imaging PRN  
Comment: Flush
- ☐ **HEParin 2000 Units in 1000 mL 0.9% NaCl infusion** Dose: 100 mL Frequency: once in imaging PRN Comment: Flush
- ☐ **sodium chloride 0.9 % bolus** Dose: 50 mL Route: intravenous Frequency: once in imaging Minimum Infusion Duration: 1.000 Minutes PRN Comment: Flush
- ☐ **sodium chloride 0.9% infusion** Dose: 250 mL Route: intravenous Frequency: once in imaging

**Medications**

**CT Contrasts**

- ☐ **barium (VOLUMEN/NEULUMEX) 0.1 % suspension** Dose: 450 mL Route: oral Frequency: once in imaging PRN  
Reasons: contrast
- ☐ **diatrizoate meglumine (HYPAQUE) 30 % solution** Dose: 300 mL Route: urethral Frequency: once in imaging PRN  
Reasons: contrast
- ☐ **iohexol (OMNIPAQUE 240) solution** Dose: 50 mL Frequency: once in imaging PRN Reasons: contrast
- ☐ **iohexol (OMNIPAQUE 300) solution** Dose: 300 Route: injection Frequency: once in imaging PRN Reasons: contrast
- ☐ **iohexol (OMNIPAQUE 300) oral solution** Dose: 30 mL Route: oral Frequency: once in imaging PRN Reasons: contrast  
Admin Instructions:  
\*\*FOR CT SCAN ONLY: Mix Iohexol 30 mL in 32 ounces of water, may add one packet of Crystal Light flavored to taste.\*\*
- ☐ **iohexol (OMNIPAQUE 350) solution** Dose: 350 Route: intravenous Frequency: once in imaging PRN Reasons: contrast
- ☐ **barium (READI-CAT 2) 2.1 % (w/v), 2.0 % (w/w) suspension** Dose: 450 mL Route: oral Frequency: once in imaging PRN  
Reasons: contrast
- ☐ **iodixanol (VISIPAQUE 320) injection** Dose: 320 Route: intravenous Frequency: once in imaging PRN Reasons: contrast
- ☐ **diatrizoate meg-diatrizoat sod (GASTROGRAFIN) 66-10 % solution** Dose: 120 mL Frequency: once in imaging PRN  
Reasons: contrast
- ☐ **iopamidol (ISOVUE-300) 61 % injection** Dose: 100 mL Route: intravenous Frequency: once in imaging PRN Reasons: contrast

**Additional Orders**

Sign: \_\_\_\_\_ Printed Name: \_\_\_\_\_ Date/Time: \_\_\_\_\_